

## FERGUSSON (POPPEN) SUCTION TUBE

11 Fr.; 45° angled

14.5cm (5.8") working length



**SKU:** 085-034-11

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Suction Tubes](#)  
**AUTHORS** [FERGUSSON, POPPEN](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)  
**TYPE** [Suction Tubes](#)  
**AUTHORS** [FERGUSSON, POPPEN](#)