

FERGUSSON (POPPEN) SUCTION TUBE

15 Fr.; 45° angled

14.5cm (5.8") working length



SKU: 085-034-15

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-----------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | FERGUSSON, POPPEN |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-----------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | FERGUSSON, POPPEN |