

FERGUSSON (POPPEN) SUCTION TUBE

5Fr.; 45° angled

14.5cm (5.8") working length



SKU: 085-034-05

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	FERGUSSON, POPPEN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	FERGUSSON, POPPEN