

## FERROZELL LIGATURE CONDUCTOR

straight

24.0cm (9.6")

**SKU:** 056-520-24

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                     |
|-------------------|-------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>                |
| <b>TYPE</b>       | <a href="#">Ligature Conductors</a> |
| <b>Aesculap</b>   | BM938                               |
| <b>KLS Martin</b> | 20-870-24                           |
| <b>AUTHORS</b>    | <a href="#">FERROZELL</a>           |

### ZUSÄTZLICHE INFORMATIONEN

|                   |                                     |
|-------------------|-------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>                |
| <b>TYPE</b>       | <a href="#">Ligature Conductors</a> |
| <b>Aesculap</b>   | BM938                               |
| <b>KLS Martin</b> | 20-870-24                           |
| <b>AUTHORS</b>    | <a href="#">FERROZELL</a>           |