

FERROZELL LIGATURE CONDUCTOR

straight

24.0cm (9.6")

SKU: 056-520-24

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM938
KLS Martin	20-870-24
AUTHORS	FERROZELL

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM938
KLS Martin	20-870-24
AUTHORS	FERROZELL