

## FIBER INTAKE/ ADAPTER

long

for Laser working elements:

630-046-08; 630-046-12; 630-047-08; 630-047-12



**SKU:** 630-047-26

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Instruments for Laser Therapy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>