

## FIBER OPTIC HEAD LAMP

focusing; Storz connection

with fiber optic cable 3.5 x 2300mm



**SKU:** 100-170-00

**Categories:** [Diagnostics](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>KLS Martin</b>	17-501-02

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Head Lamps</a>
<b>KLS Martin</b>	17-501-02