

## FIBER OPTIC LIGHT CARRIER

for operating laryngoscopes fig.3-6



**SKU:** 103-404-00

**Categories:** [Anaesthesia](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Operation Laryngoscopes](#)

**Aesculap** OP297

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Operation Laryngoscopes](#)

**Aesculap** OP297