

## FIBER OPTIC LIGHT CARRIER

for operating laryngoscopes fig.1-2

**SKU:** 103-403-00

**Categories:** [Anaesthesia](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Operation Laryngoscopes</a>
<b>Aesculap</b>	OP298

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Operation Laryngoscopes</a>
<b>Aesculap</b>	OP298