

## FIBRINADHESIVE NEEDLE

tube 2.6mm, 230cm, w/ mixing chamber

needle, 0.9 x 5mm (single use)



**SKU:** PK-8KZEW2609230-S

**Category:** [Endoscopy](#)



## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

## ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)