

FINOCHIETTO NEEDLE HOLDER

angled; serrated; grooved jaw

27.0cm (10.8")




SKU: 050-211-27

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM277
KLS Martin	20-415-27
AUTHORS	FINOCHIETTO
Medicon	10.43.27
SPECIFICATIONS	angled , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM277
KLS Martin	20-415-27

AUTHORS

[FINOCHIETTO](#)

Medicon

10.43.27

SPECIFICATIONS

[angled](#), [serrated](#)