

FINSEN SELFRETAINING RETRACTOR

sharp; 3X4 teeth

7.0cm (2.8")




SKU: 070-610-07

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV002
KLS Martin	15-702-07
AUTHORS	FINSEN
Medicon	24.04.07
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV002
KLS Martin	15-702-07

AUTHORS

[FINSEN](#)

Medicon

24.04.07

SPECIFICATIONS

[sharp, toothed](#)