

FINSEN SELFRETAINING RETRACTOR

blunt; 2X3 teeth

5.0cm (2")


SKU: 070-611-05

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-701-05
AUTHORS	FINSEN
Medicon	24.04.06
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-701-05
AUTHORS	FINSEN
Medicon	24.04.06

SPECIFICATIONS

[blunt](#), [toothed](#)