

FINSEN SELFRETAINING RETRACTOR

blunt; 3X4 teeth

7.0cm (2.8")




SKU: 070-611-07

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV003
KLS Martin	15-703-07
AUTHORS	FINSEN
Medicon	24.04.08
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV003
KLS Martin	15-703-07

AUTHORS

[FINSEN](#)

Medicon

24.04.08

SPECIFICATIONS

[blunt](#), [toothed](#)