

FINSEN SELFRETAINING RETRACTOR

sharp; 2X3 teeth

5.0cm (2")


SKU: 070-610-05

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV001
KLS Martin	15-700-05
AUTHORS	FINSEN
Medicon	24.04.05
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV001
KLS Martin	15-700-05

AUTHORS

[FINSEN](#)

Medicon

24.04.05

SPECIFICATIONS

[sharp, toothed](#)