

## FISCHER TENDON STRIPPER

Ø 4.0mm; straight

31.0cm (12.4")

**SKU:** 195-262-31

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>KLS Martin</b>	23-660-31
<b>AUTHORS</b>	<a href="#">FISCHER</a>
<b>Medicon</b>	78.22.30
<b>SPECIFICATIONS</b>	<a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tendon Stripper</a>
<b>KLS Martin</b>	23-660-31
<b>AUTHORS</b>	<a href="#">FISCHER</a>
<b>Medicon</b>	78.22.30

## SPECIFICATIONS

[straight](#)