

## FLUTED BIT ELECTRODE; BIPOAR; INSERT

bit 1mm; for Ø 5.0mm shaft

for 628-201-53 and 628-205-53

**SKU:** 628-242-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>