

## FLUTED BIT ELECTRODE BIPOLAR; 1MM; INSERT

short insulated; for coagulation forceps

628-101-00 and 628-105-00; 330mm

No CE mark! No labelling!



**SKU:** 628-142-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">HF-Instruments, bipolar</a>
<b>TYPE</b>	<a href="#">Electrodes</a>