

## FOERSTER-BALLENGER SPONGE FORCEPS

straight; smooth

18.0cm (7.2")




**SKU:** 030-170-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>AUTHORS</b>	<a href="#">BALLENGER</a> , <a href="#">FOERSTER</a>
<b>Medicon</b>	19.52.18
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Sponge and Dressing Forceps</a>
<b>AUTHORS</b>	<a href="#">BALLENGER</a> , <a href="#">FOERSTER</a>
<b>Medicon</b>	19.52.18
<b>SPECIFICATIONS</b>	<a href="#">straight</a>