

## FOMON ALAR HOOKLET

2-pronged; button ended

16.0cm (6.4")


**SKU:** 112-305-16

**Categories:** [Open Surgery](#), [Rhinology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	OL202
<b>KLS Martin</b>	37-706-16
<b>AUTHORS</b>	<a href="#">FOMON</a>
<b>SPECIFICATIONS</b>	<a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	OL202
<b>KLS Martin</b>	37-706-16
<b>AUTHORS</b>	<a href="#">FOMON</a>

## SPECIFICATIONS

[toothed](#)