

FOREIGN BODY FORCEPS; Ø 1.7MM; FL 600MM

spoon serrated; long jaw part; semi-rigid

double action; insert only



SKU: 583-472-16A

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments