

## FORIGN BODY FORCEPS; Ø 1.7MM; FL 600MM

spoon serrated; long jaw part;semi-rigid

double action; with LuerLock;insert only



with Luer Lock

**SKU:** 583-472-16LLA

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Uretero-Renoscapy](#)

**TYPE** [Semi-rigid Instruments](#)