

FRAZIER SUCTION TUBE

15 Fr.; 45° angled

13.0cm (5.2") working length

SKU: 085-031-15

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	FRAZIER
Medicon	09.83.34

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	FRAZIER
Medicon	09.83.34