

FRAZIER SUCTION TUBE

12 Fr.; 45° angled

18.0cm (7.2") working length

SKU: 085-033-12

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF374
KLS Martin	24-769-17
AUTHORS	FRAZIER
Medicon	09.83.82

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF374
KLS Martin	24-769-17
AUTHORS	FRAZIER

Medicon

09.83.82