

FRAZIER SUCTION TUBE

10 Fr.; 45° angled

10.0cm (4") working length

SKU: 085-030-10

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF923
KLS Martin	18-523-10
AUTHORS	FRAZIER
Medicon	09.83.50

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF923
KLS Martin	18-523-10
AUTHORS	FRAZIER

Medicon

09.83.50