

## FUCHSIG HAEMOSTATIC FORCEPES

curved; serrated

20.0cm (8")



**SKU:** 040-419-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH198
<b>KLS Martin</b>	13-369-20
<b>AUTHORS</b>	<a href="#">FUCHSIG</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH198
<b>KLS Martin</b>	13-369-20
<b>AUTHORS</b>	<a href="#">FUCHSIG</a>

## SPECIFICATIONS

[curved](#), [serrated](#)