

GELPI-LOKTITE SELFRETAINING RETRACTOR

sharp

18.0cm (7.2")




SKU: 070-250-19

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-804-18
AUTHORS	GELPI , LOKTITE
Medicon	24.40.57
SPECIFICATIONS	sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
KLS Martin	15-804-18
AUTHORS	GELPI , LOKTITE
Medicon	24.40.57

SPECIFICATIONS

[sharp](#)