

## GELPI-LOKTITE SELFRETAINING RETRACTOR

sharp

18.0cm (7.2")




**SKU:** 070-250-19

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-804-18
<b>AUTHORS</b>	<a href="#">GELPI</a> , <a href="#">LOKTITE</a>
<b>Medicon</b>	24.40.57
<b>SPECIFICATIONS</b>	<a href="#">sharp</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>KLS Martin</b>	15-804-18
<b>AUTHORS</b>	<a href="#">GELPI</a> , <a href="#">LOKTITE</a>
<b>Medicon</b>	24.40.57

## SPECIFICATIONS

[sharp](#)