

GIL-VERNET RETRACTOR

11mm; flexible

24.0cm (9.6")

SKU: 062-230-11

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-11
AUTHORS	GIL-VERNET



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-11
AUTHORS	GIL-VERNET