

GIL-VERNET RETRACTOR

20mm; flexible

24.0cm (9.6")


SKU: 062-230-20

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-20
AUTHORS	GIL-VERNET



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-20
AUTHORS	GIL-VERNET