

## GIL-VERNET RETRACTOR

23mm; flexible

24.0cm (9.6")


**SKU:** 062-230-23

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-174-23
<b>AUTHORS</b>	<a href="#">GIL-VERNET</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-174-23
<b>AUTHORS</b>	<a href="#">GIL-VERNET</a>