

GIL-VERNET RETRACTOR

8mm; flexible

24.0cm (9.6")

SKU: 062-230-08

Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Delicate Hooks |
| KLS Martin | 15-174-08 |
| AUTHORS | GIL-VERNET |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Delicate Hooks |
| KLS Martin | 15-174-08 |
| AUTHORS | GIL-VERNET |