

GLOVER BULLDOG CLAMP

straight; atraumatic

7.0cm (2.8")

SKU: 026-298-06

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB432
KLS Martin	24-398-06
AUTHORS	GLOVER
SPECIFICATIONS	atraumatic , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB432
KLS Martin	24-398-06
AUTHORS	GLOVER

SPECIFICATIONS

[atraumatic](#), [straight](#)