

GLOVER BULLDOG CLAMP

straight; atraumatic

8.0cm (3.2")

SKU: 026-298-07

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB432
KLS Martin	24-398-06
AUTHORS	GLOVER
Medicon	15.18.07
SPECIFICATIONS	atraumatic , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB432
KLS Martin	24-398-06

AUTHORS

[GLOVER](#)

Medicon

15.18.07

SPECIFICATIONS

[atraumatic](#), [straight](#)