

## GLOVER BULLDOG CLAMP

straight; atraumatic

8.0cm (3.2")

**SKU:** 026-298-07

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB432
<b>KLS Martin</b>	24-398-06
<b>AUTHORS</b>	<a href="#">GLOVER</a>
<b>Medicon</b>	15.18.07
<b>SPECIFICATIONS</b>	<a href="#">atraumatic</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB432
<b>KLS Martin</b>	24-398-06

## AUTHORS

[GLOVER](#)

Medicon

15.18.07

## SPECIFICATIONS

[atraumatic](#), [straight](#)