

GLOVER BULLDOG CLAMP

straight; atraumatic

11.0cm (4.4")

SKU: 026-298-11

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB436
KLS Martin	24-398-11
AUTHORS	GLOVER
Medicon	15.18.11
SPECIFICATIONS	atraumatic , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB436
KLS Martin	24-398-11

AUTHORS

[GLOVER](#)

Medicon

15.18.11

SPECIFICATIONS

[atraumatic](#), [straight](#)