

## GLOVER BULLDOG CLAMP

angled; atraumatic

6.5cm (2.6")



**SKU:** 026-193-06

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB431
<b>KLS Martin</b>	24-293-06
<b>AUTHORS</b>	<a href="#">GLOVER</a>
<b>Medicon</b>	15.19.55
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">atraumatic</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	FB431
<b>KLS Martin</b>	24-293-06

## AUTHORS

[GLOVER](#)

## Medicon

15.19.55

## SPECIFICATIONS

[angled](#), [atraumatic](#)