

GRAEFE IRIS FORCEPS

straight; serrated

7.0cm (2.8")



SKU: 020-810-07

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Iris Forceps
KLS Martin	35-360-07
AUTHORS	GRAEFE
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Iris Forceps
KLS Martin	35-360-07
AUTHORS	GRAEFE
SPECIFICATIONS	serrated , straight