

## GRASPING FORCEP ELECTRODE;BIPOLAR;INSERT

fenestrated; for Ø 5.0mm shaft

for article 628-301-53 and 628-305-53

**SKU:** 628-344-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                    |
| <b>RANGE</b> | <a href="#">HF-Instruments, bipolar</a> |
| <b>TYPE</b>  | <a href="#">Electrodes</a>              |

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                    |
| <b>RANGE</b> | <a href="#">HF-Instruments, bipolar</a> |
| <b>TYPE</b>  | <a href="#">Electrodes</a>              |