

## GRASPING FORCEPS, FENESTRATED ELECTRODE

for bipolar coagulation forceps

290mm



**SKU:** 628-444-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Pediatric Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Pediatric Laparoscopy](#)