

## GRAY GALL DUCT CLAMP

s-shaped; curved; serrated

22.0cm (8.8")




**SKU:** 040-523-22

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Gall Duct Forceps</a>
<b>Aesculap</b>	EB011
<b>KLS Martin</b>	13-561-22
<b>AUTHORS</b>	<a href="#">GRAY</a>
<b>Medicon</b>	51.37.01
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Gall Duct Forceps</a>
<b>Aesculap</b>	EB011
<b>KLS Martin</b>	13-561-22

## AUTHORS

[GRAY](#)

## Medicon

51.37.01

## SPECIFICATIONS

[curved](#), [serrated](#)