

GREGORY-SOFT CAROTID CLAMP

angled; atraumatic

11.0cm (4.4")

SKU: 026-197-11

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	24-297-11
AUTHORS	GREGORY , SOFT
SPECIFICATIONS	angled , atraumatic



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	24-297-11
AUTHORS	GREGORY , SOFT
SPECIFICATIONS	angled , atraumatic