

## GREGORY-SOFT CAROTID CLAMP

angled; atraumatic

9.0cm (3.6")

**SKU:** 026-197-09

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	24-297-09
<b>AUTHORS</b>	<a href="#">GREGORY</a> , <a href="#">SOFT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">atraumatic</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	24-297-09
<b>AUTHORS</b>	<a href="#">GREGORY</a> , <a href="#">SOFT</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">atraumatic</a>