

## GROOVED DIRECTOR

Ø 3.5mm

20.0cm (8")

**SKU:** 080-110-20

**Categories:** [Open Surgery](#), [Probes](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Grooved Director</a>
<b>KLS Martin</b>	16-160-20
<b>Medicon</b>	09.15.20

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Grooved Director</a>
<b>KLS Martin</b>	16-160-20
<b>Medicon</b>	09.15.20