

GRUENWALD NASAL TAMPON FORCEPS

bayonet shaped; serrated

20.0cm (8"); width of jaw 2.6 +/-0.2



SKU: 020-790-20

Categories: [Forceps](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing and Tampon Forceps, Cover Glass Forceps
Aesculap	BD883
KLS Martin	37-115-20
AUTHORS	GRAEFE
Medicon	07.30.19
SPECIFICATIONS	bajonett , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing and Tampon Forceps, Cover Glass Forceps
Aesculap	BD883
KLS Martin	37-115-20

AUTHORS

[GRAEFE](#)

Medicon

07.30.19

SPECIFICATIONS

[bajonett](#), [serrated](#)