

## GUILLEN SUCTION RASPATORY

right curved tip; Ø 4.0mm; 20.0cm

**SKU:** 085-076-20

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">GUILLEN</a>       |

### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">GUILLEN</a>       |