

## GUILLEN SUCTION RASPATORY

straight tip; Ø 4.0mm; 20.0cm

**SKU:** 085-074-20

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	37-840-20
<b>AUTHORS</b>	<a href="#">GUILLEN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	37-840-20
<b>AUTHORS</b>	<a href="#">GUILLEN</a>