

GUNNAR HEY TUBING CLAMP

Roller clamp

18.5cm (7.4")



SKU: 032-090-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	GUNNAR HEY
SPECIFICATIONS	straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
AUTHORS	GUNNAR HEY
SPECIFICATIONS	straight