

## GUNNAR HEY TUBING CLAMP

Roller clamp

18.5cm (7.4")



**SKU:** 032-090-18

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">GUNNAR HEY</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">GUNNAR HEY</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>