

GUYON KIDNEY PEDICLE CLAMP

curved; serrated

22.0cm (8.8")

SKU: 040-509-22

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH837
KLS Martin	13-531-22
AUTHORS	GUYON
Medicon	51.56.22
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH837
KLS Martin	13-531-22

AUTHORS

[GUYON](#)

Medicon

51.56.22

SPECIFICATIONS

[curved](#), [serrated](#)